

Product Datasheet

AKAP9 Antibody (orb1242533)

Description

AKAP9 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human AKAP9.

Antibody used in WB on Human Jurkat 1.25...

Target

AKAP9

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store AKAP9 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

AKAP9 antibody can be used for detection of AKAP9 by ELISA at 1:12500. AKAP9 antibody can be used for detection of AKAP9 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

36 kDa, 453 kDa, 452 kDa

Uniprot ID

[Q6PJH3](#)

NCBI

[AAH15533](#)

Dilution Range

AKAP9 antibody can be used for detection of AKAP9 by ELISA at 1:12500. AKAP9 antibody can be used for detection of AKAP9 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.