


## Product Datasheet

### TAF1C Antibody (orb1242536)

| Description                | TAF1C Antibody                                                                                                                                                                                                                     |                                                                                     |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                             |  |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                                              |                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                       |                                                                                     |
| <b>Tested Applications</b> | ELISA, WB                                                                                                                                                                                                                          |                                                                                     |
| <b>Immunogen</b>           | Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TAF1C.                                                                                                                             |                                                                                     |
| <b>Target</b>              | TAF1C                                                                                                                                                                                                                              | Antibody used in WB on Human Jurkat cell...                                         |
| <b>Preservatives</b>       | Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.                                                                                                                                          |                                                                                     |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                                                                             |                                                                                     |
| <b>Concentration</b>       | batch dependent                                                                                                                                                                                                                    |                                                                                     |
| <b>Storage</b>             | For short periods of storage (days) store at 4°C. For longer periods of storage, store TAF1C antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.                                                              |                                                                                     |
| <b>Note</b>                | For research use only                                                                                                                                                                                                              |                                                                                     |
| <b>Application notes</b>   | TAF1C antibody can be used for detection of TAF1C by ELISA at 1:12500. TAF1C antibody can be used for detection of TAF1C by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000. |                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                         |                                                                                     |
| <b>MW</b>                  | 85 kDa, 95 kDa                                                                                                                                                                                                                     |                                                                                     |
| <b>Uniprot ID</b>          | <a href="#">Q15572</a>                                                                                                                                                                                                             |                                                                                     |
| <b>NCBI</b>                | <a href="#">NP_647610</a>                                                                                                                                                                                                          |                                                                                     |
| <b>Dilution Range</b>      | TAF1C antibody can be used for detection of TAF1C by ELISA at 1:12500. TAF1C antibody can be used for detection of TAF1C by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000. |                                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                    |                                                                                     |