

Product Datasheet

CIR1 Antibody (orb1242566)

Description

CIR1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CIR.

Target

CIR1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CIR antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CIR antibody can be used for detection of CIR by ELISA at 1:62500. CIR antibody can be used for detection of CIR by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

52 kDa

Uniprot ID

[Q86X95](#)

NCBI


[NP_004873](#)

Dilution Range

CIR antibody can be used for detection of CIR by ELISA at 1:62500. CIR antibody can be used for detection of CIR by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody
used in WB
on
Transfected
293T ...