

Product Datasheet

ETV5 Antibody (orb1242576)

Description

ETV5 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ETV5.

Target

ETV5

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ETV5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ETV5 antibody can be used for detection of ETV5 by ELISA at 1:62500. ETV5 antibody can be used for detection of ETV5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

58 kDa

Uniprot ID

[P41161](#)

NCBI

[NP_004445](#)

Dilution Range

ETV5 antibody can be used for detection of ETV5 by ELISA at 1:62500. ETV5 antibody can be used for detection of ETV5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

81kDa

40kDa

31kDa

22kDa

Antibody used in WB on Transfected 293T ...

ZBTB48

90

65

40

31

22

Rabbit Anti-ZBTB48
Sample Type: Human 293T
Antibody Concentration: 1µg/mL

Antibody used in WB on Hum. Fetal Brain ...

ETV5

90

65

40

31

22

Rabbit Anti-ETV5
Sample Type: Human Fetal Heart
Antibody Concentration: 1µg/mL

Antibody used in WB on Hum. Fetal Heart ...