

Product Datasheet

RNF14 Antibody (orb1242577)

Description

RNF14 Antibody



Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Antibody used in IHC on Human kidney.

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RNF14.

Target

RNF14

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store RNF14 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in WB on Human Jurkat 2.5 ...

Note

For research use only

Application notes

RNF14 antibody can be used for detection of RNF14 by ELISA at 1:312500. RNF14 antibody can be used for detection of RNF14 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

52 kDa, 54 kDa, 40 kDa, 54 kDa, 54 kDa

Uniprot ID

[Q9UBS8](#)

NCBI

[NP_899648](#)

Dilution Range

RNF14 antibody can be used for detection of RNF14 by ELISA at 1:312500. RNF14 antibody can be used for detection of RNF14 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.