

Product Datasheet

TIMELESS Antibody (orb1242586)

Description

TIMELESS Antibody

Species/Host	Rabbit
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TIMELESS.
Target	TIMELESS
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store TIMELESS antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	TIMELESS antibody can be used for detection of TIMELESS by ELISA at 1:312500. TIMELESS antibody can be used for detection of TIMELESS by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	139 kDa
Uniprot ID	Q9UNS1
NCBI	NP_003911
Dilution Range	TIMELESS antibody can be used for detection of TIMELESS by ELISA at 1:312500. TIMELESS antibody can be used for detection of TIMELESS by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

Antibody
used in WB
on
Transfected
293T ...