

Product Datasheet


PSMC5 Antibody (orb1242631)

Description

PSMC5 Antibody


Species/Host	Rabbit
Reactivity	C. elegans, Canine, Drosophila, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PSMC5.
Target	PSMC5
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PSMC5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	PSMC5 antibody can be used for detection of PSMC5 by ELISA at 1:312500. PSMC5 antibody can be used for detection of PSMC5 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	46 kDa
Uniprot ID	P62195
NCBI	NP_002796
Dilution Range	PSMC5 antibody can be used for detection of PSMC5 by ELISA at 1:312500. PSMC5 antibody can be used for detection of PSMC5 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

PSMC5



Antibody used in WB on Mouse Brain at 2 ...

PSMC5



See Other Application 1 Data and Customer Feedback for more information

Antibody used in IP on Mouse brain at 1:...



Antibody used in WB on Human Jurkat 1.25...