

## Product Datasheet

### PBX2 Antibody (orb1242634)

## Description

PBX2 Antibody

### Species/Host

Rabbit

### Reactivity

Canine, Human, Mouse, Rat, Zebrafish

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

### Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PBX2.

### Target

PBX2

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PBX2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

### Note

For research use only

### Application notes

PBX2 antibody can be used for detection of PBX2 by ELISA at 1:62500. PBX2 antibody can be used for detection of PBX2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

46 kDa

### Uniprot ID

[P40425](#)

### NCBI

[NP\\_002577](#)

### Dilution Range

PBX2 antibody can be used for detection of PBX2 by ELISA at 1:62500. PBX2 antibody can be used for detection of PBX2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.



Antibody  
used in WB  
on Human  
Jurkat 1.25...