

Product Datasheet

NR1H5 Antibody (orb1242735)

Description

NR1H5 Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Mouse
SP2/0 at 2....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of mouse NR1H5.

Target

NR1H5

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store NR1H5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

NR1H5 antibody can be used for detection of NR1H5 by ELISA at 1:312500. NR1H5 antibody can be used for detection of NR1H5 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

49 kDa

Uniprot ID

[Q811W9](#)

NCBI

[NP_963916](#)

Dilution Range

NR1H5 antibody can be used for detection of NR1H5 by ELISA at 1:312500. NR1H5 antibody can be used for detection of NR1H5 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.