

Product Datasheet

Duxbl1 Antibody (orb1242742)

Description		
Species/Host	Rabbit	
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of mouse 1110051B16RIK.	Antibody used in WB on Mouse SP2/0 at 0....
Target	Duxbl1	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store 1110051B16RIK antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	1110051B16RIK antibody can be used for detection of 1110051B16RIK by ELISA at 1:1562500. 1110051B16RIK antibody can be used for detection of 1110051B16RIK by western blot at 0.125 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	39 kDa	
Uniprot ID	Q7TNE6	
NCBI	NP_899245	
Dilution Range	1110051B16RIK antibody can be used for detection of 1110051B16RIK by ELISA at 1:1562500. 1110051B16RIK antibody can be used for detection of 1110051B16RIK by western blot at 0.125 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	