

Product Datasheet

Radil Antibody (orb1242746)

Description	Radil Antibody	
Species/Host	Rabbit	
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Mouse SP2/0 at 5....
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of mouse D930005D10RIK.	
Target	Radil	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store D930005D10RIK antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	D930005D10RIK antibody can be used for detection of D930005D10RIK by ELISA at 1:62500. D930005D10RIK antibody can be used for detection of D930005D10RIK by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	109 kDa	
Uniprot ID	Q69Z89	
NCBI	NP_848817	
Dilution Range	D930005D10RIK antibody can be used for detection of D930005D10RIK by ELISA at 1:62500. D930005D10RIK antibody can be used for detection of D930005D10RIK by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	