


Product Datasheet

E2f7 Antibody (orb1242748)


Description

E2f7 Antibody

Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of mouse E2f7.
Target	E2f7
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store E2f7 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	E2f7 antibody can be used for detection of E2f7 by ELISA at 1:62500. E2f7 antibody can be used for detection of E2f7 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	99 kDa
Uniprot ID	Q8BSQ3
NCBI	NP_848724
Dilution Range	E2f7 antibody can be used for detection of E2f7 by ELISA at 1:62500. E2f7 antibody can be used for detection of E2f7 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Mouse stomach.



Antibody used in WB on Mouse SP2/0 at 0....