

Product Datasheet

Cphx1 Antibody (orb1242752)

Description

Cphx1 Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Mouse
SP2/0 at 5....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of mouse C330003B14RIK.

Target

Cphx1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store C330003B14RIK antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

C330003B14RIK antibody can be used for detection of C330003B14RIK by ELISA at 1:62500. C330003B14RIK antibody can be used for detection of C330003B14RIK by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

20 kDa

Uniprot ID

[Q8BX39](#)

NCBI

[NP_780551](#)

Dilution Range

C330003B14RIK antibody can be used for detection of C330003B14RIK by ELISA at 1:62500. C330003B14RIK antibody can be used for detection of C330003B14RIK by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.