

Product Datasheet

OBOX6 Antibody (orb1242761)

Description		
Species/Host	Rabbit	
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of mouse OBOX6.	Antibody used in WB on Mouse SP2/0 at 1....
Target	OBOX6	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store OBOX6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	OBOX6 antibody can be used for detection of OBOX6 by ELISA at 1:1562500. OBOX6 antibody can be used for detection of OBOX6 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	34 kDa	
Uniprot ID	Q8VHG3	
NCBI	NP_663756	
Dilution Range	OBOX6 antibody can be used for detection of OBOX6 by ELISA at 1:1562500. OBOX6 antibody can be used for detection of OBOX6 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	