

Product Datasheet

Ces2a Antibody (orb1242765)

Description

Ces2a Antibody

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Mouse
NIH-3T3 at ...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of mouse CES6.

Target

Ces2a

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CES6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CES6 antibody can be used for detection of CES6 by ELISA at 1:312500. CES6 antibody can be used for detection of CES6 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

61 kDa

Uniprot ID

[Q8QZR3](#)

NCBI

[NP_598721](#)

Dilution Range

CES6 antibody can be used for detection of CES6 by ELISA at 1:312500. CES6 antibody can be used for detection of CES6 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.