

Product Datasheet


NFAT5 Antibody (orb1242766)

Description

NFAT5 Antibody

Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of mouse NFAT5.
Target	NFAT5
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store NFAT5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	NFAT5 antibody can be used for detection of NFAT5 by ELISA at 1:62500. NFAT5 antibody can be used for detection of NFAT5 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	74 kDa
Uniprot ID	Q9JIH4
NCBI	AAF79117
Dilution Range	NFAT5 antibody can be used for detection of NFAT5 by ELISA at 1:62500. NFAT5 antibody can be used for detection of NFAT5 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

NFAT5
1 2 3



See datasheet 12000 and Customer Feedback for more information

Antibody used in WB on HT-29 at 1:1000 (...)



Antibody used in WB on Mouse SP2/0 at 0....