

Product Datasheet

Tsc22d4 Antibody (orb1242792)

Description

Tsc22d4 Antibody

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of mouse Tsc22d4.

Target

Tsc22d4

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store Tsc22d4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

Tsc22d4 antibody can be used for detection of Tsc22d4 by ELISA at 1:62500. Tsc22d4 antibody can be used for detection of Tsc22d4 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

43 kDa

Uniprot ID

[Q9EQN3](#)

NCBI


[NP_076399](#)

Dilution Range

Tsc22d4 antibody can be used for detection of Tsc22d4 by ELISA at 1:62500. Tsc22d4 antibody can be used for detection of Tsc22d4 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Mouse stomach.



Antibody used in WB on Mouse SP2/0 at 5....