

Product Datasheet

Nkx3-2 Antibody (orb1242874)

Description

Nkx3-2 Antibody

Species/Host

Rabbit

Reactivity

Canine, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of mouse Bapx1.

Target

Nkx3-2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store Bapx1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

Bapx1 antibody can be used for detection of Bapx1 by ELISA at 1:62500. Bapx1 antibody can be used for detection of Bapx1 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

37 kDa

Uniprot ID

[P97503](#)

NCBI

[NP_031550](#)

Dilution Range

Bapx1 antibody can be used for detection of Bapx1 by ELISA at 1:62500. Bapx1 antibody can be used for detection of Bapx1 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.


70 kDa

40 kDa

35 kDa

21 kDa

Antibody used in WB on Human Jurkat 0.2-...



Antibody used in WB on Jurkat 1 µg/ml (L...