

Product Datasheet

KCTD11 Antibody (orb1242878)

Description

KCTD11 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KCTD11.

Target

KCTD11

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store KCTD11 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

KCTD11 antibody can be used for detection of KCTD11 by ELISA at 1:62500. KCTD11 antibody can be used for detection of KCTD11 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

26 kDa

Uniprot ID

[Q693B1](#)

NCBI

[NP_001002914](#)

Dilution Range

KCTD11 antibody can be used for detection of KCTD11 by ELISA at 1:62500. KCTD11 antibody can be used for detection of KCTD11 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

98 kDa

42 kDa

32 kDa

27 kDa

Antibody used in WB on Human HepG2 at 1....

KCTD11

98

65

40

31

22

Rabbit Anti-KCTD11
Sample Type: 721_B
Antibody Concentration: 1 µg/mL

Antibody used in WB on Human 721_B at 1 ...