

Product Datasheet

WRNIP1 Antibody (orb1242944)

Description

WRNIP1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human WRNIP1.

Target

WRNIP1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store WRNIP1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

WRNIP1 antibody can be used for detection of WRNIP1 by ELISA at 1:62500. WRNIP1 antibody can be used for detection of WRNIP1 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

72 kDa

Uniprot ID

[Q96S55](#)

NCBI

[NP_064520](#)

Dilution Range

WRNIP1 antibody can be used for detection of WRNIP1 by ELISA at 1:62500. WRNIP1 antibody can be used for detection of WRNIP1 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

87kDa

66kDa


46kDa

36kDa

21kDa

A B

Antibody used in WB on Human HepG2 at 0....



Antibody used in WB on Human at 1:1000.