

Product Datasheet

DDX49 Antibody (orb1242945)

Description

DDX49 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Thymus at
2...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human DDX49.

Target

DDX49

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store DDX49 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

DDX49 antibody can be used for detection of DDX49 by ELISA at 1:312500. DDX49 antibody can be used for detection of DDX49 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

53 kDa

Uniprot ID

[Q9Y6V7](#)

NCBI

[NP_061943](#)

Dilution Range

DDX49 antibody can be used for detection of DDX49 by ELISA at 1:312500. DDX49 antibody can be used for detection of DDX49 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.