

## Product Datasheet

### TGIF2LY Antibody (orb1243001)

## Description

## TGIF2LY Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

### Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TGIF2LY.

### Target

TGIF2LY

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TGIF2LY antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

### Note

For research use only

### Application notes

TGIF2LY antibody can be used for detection of TGIF2LY by ELISA at 1:62500. TGIF2LY antibody can be used for detection of TGIF2LY by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

21 kDa, 27 kDa, 27 kDa

### Uniprot ID

[Q8IUE0](#)

### NCBI

[NP\\_631960](#)

### Dilution Range

TGIF2LY antibody can be used for detection of TGIF2LY by ELISA at 1:62500. TGIF2LY antibody can be used for detection of TGIF2LY by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.



Antibody used in WB on Human HepG2 at 1....