

Product Datasheet

ZNF551 Antibody (orb1243003)

Description

ZNF551 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZNF551.

Target

ZNF551

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZNF551 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ZNF551 antibody can be used for detection of ZNF551 by ELISA at 1:1562500. ZNF551 antibody can be used for detection of ZNF551 by western blot at 0.6 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

40 kDa, 76 kDa

Uniprot ID

[Q7Z340](#)

NCBI



[BAC03625](#)

Dilution Range


ZNF551 antibody can be used for detection of ZNF551 by ELISA at 1:1562500. ZNF551 antibody can be used for detection of ZNF551 by western blot at 0.6 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Liver.



Antibody used in WB on Human Jurkat 0.6 ...