


Product Datasheet

HES4 Antibody (orb1243037)

Description	HES4 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Transfected 293T ...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human HES4.	
Target	HES4	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store HES4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	HES4 antibody can be used for detection of HES4 by ELISA at 1:312500. HES4 antibody can be used for detection of HES4 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	24 kDa	
Uniprot ID	Q9HCC6	
NCBI	NP_066993	
Dilution Range	HES4 antibody can be used for detection of HES4 by ELISA at 1:312500. HES4 antibody can be used for detection of HES4 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	