

Product Datasheet

ZNF223 Antibody (orb1243067)

Description

ZNF223 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZNF223.

Target

ZNF223

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZNF223 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ZNF223 antibody can be used for detection of ZNF223 by ELISA at 1:12500. ZNF223 antibody can be used for detection of ZNF223 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

56 kDa

Uniprot ID

[Q9UK11](#)

NCBI


[NP_037493](#)

Dilution Range

ZNF223 antibody can be used for detection of ZNF223 by ELISA at 1:12500. ZNF223 antibody can be used for detection of ZNF223 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Skin.



Antibody used in WB on Human Jurkat 0.2-...