

Product Datasheet

ZNF79 Antibody (orb1243073)

Description	ZNF79 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human Jurkat 0.25...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZNF79.	
Target	ZNF79	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ZNF79 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	ZNF79 antibody can be used for detection of ZNF79 by ELISA at 1:62500. ZNF79 antibody can be used for detection of ZNF79 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	55 kDa	
Uniprot ID	Q15937	
NCBI	NP_009066	
Dilution Range	ZNF79 antibody can be used for detection of ZNF79 by ELISA at 1:62500. ZNF79 antibody can be used for detection of ZNF79 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	