


Product Datasheet

ILF2 Antibody (orb1243102)

Description

ILF2 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ILF2.
Target	ILF2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ILF2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	ILF2 antibody can be used for detection of ILF2 by ELISA at 1:62500. ILF2 antibody can be used for detection of ILF2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	Q12905
NCBI	NP_004506
Dilution Range	ILF2 antibody can be used for detection of ILF2 by ELISA at 1:62500. ILF2 antibody can be used for detection of ILF2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in WB on Human Jurkat 1.25...



Antibody used in WB on Human MCF-7 at 1....