

Product Datasheet

KCNK13 Antibody (orb1243162)

Description

KCNK13 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KCNK13.

Target

KCNK13

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store KCNK13 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

KCNK13 antibody can be used for detection of KCNK13 by ELISA at 1:1562500. KCNK13 antibody can be used for detection of KCNK13 by western blot at 0.12 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

45 kDa

Uniprot ID

[Q9HB14](#)

NCBI



[NP_071337](#)

Dilution Range

KCNK13 antibody can be used for detection of KCNK13 by ELISA at 1:1562500. KCNK13 antibody can be used for detection of KCNK13 by western blot at 0.12 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.


Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Brain at 4...

KCNK13



Blue/DAPI Green/KCNK13

[See IHC 2 Data and Customer Feedback for more information](#)

Antibody used in IHC.



Antibody used in WB on Human Jurkat 0.12...