

Product Datasheet

CLIC4 Antibody (orb1243183)

Description

CLIC4 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Jurkat 2.5 ...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CLIC4 Channel4.

Target

CLIC4

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CLIC4 Channel4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CLIC4 Channel4 antibody can be used for detection of CLIC4 Channel4 by ELISA at 1:62500. CLIC4 Channel4 antibody can be used for detection of CLIC4 Channel4 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

28 kDa

Uniprot ID

[Q5VSX5](#)

NCBI

[NP_039234](#)

Dilution Range

CLIC4 Channel4 antibody can be used for detection of CLIC4 Channel4 by ELISA at 1:62500. CLIC4 Channel4 antibody can be used for detection of CLIC4 Channel4 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.