

Product Datasheet


P2RX1 Antibody (orb1243195)

Description

P2RX1 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human P2RX1.
Target	P2RX1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store P2RX1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	P2RX1 antibody can be used for detection of P2RX1 by ELISA at 1:312500. P2RX1 antibody can be used for detection of P2RX1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	45 kDa
Uniprot ID	P51575
NCBI	NP_002549
Dilution Range	P2RX1 antibody can be used for detection of P2RX1 by ELISA at 1:312500. P2RX1 antibody can be used for detection of P2RX1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

Antibody used in WB on Human HepG2 at 2....



Antibody used in IHC on Human Pancreas.