

Product Datasheet


KCNN2 Antibody (orb1243198)

Description

KCNN2 Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KCNN2.
Target	KCNN2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store KCNN2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	KCNN2 antibody can be used for detection of KCNN2 by ELISA at 1:62500. KCNN2 antibody can be used for detection of KCNN2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	64 kDa
Uniprot ID	Q9H2S1
NCBI	NP_067627
Dilution Range	KCNN2 antibody can be used for detection of KCNN2 by ELISA at 1:62500. KCNN2 antibody can be used for detection of KCNN2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

KCNN2



See HC2 Data and Customer Feedback for more information

Antibody used in IHC on Rat brain section...



Green: KCNN2
See HC1 Data and Customer Feedback for more information

Antibody used in IHC on Rhesus macaque s...



Antibody used in WB on Human HepG2 at 1....