

Product Datasheet

CLIC1 Antibody (orb1243203)

Description

CLIC1 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human CLIC1 Channel1.

Target

CLIC1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CLIC1 Channel1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CLIC1 Channel1 antibody can be used for detection of CLIC1 Channel1 by ELISA at 1:312500. CLIC1 Channel1 antibody can be used for detection of CLIC1 Channel1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

27 kDa

Uniprot ID

[Q502X1](#)

NCBI


[NP_001279](#)

Dilution Range

CLIC1 Channel1 antibody can be used for detection of CLIC1 Channel1 by ELISA at 1:312500. CLIC1 Channel1 antibody can be used for detection of CLIC1 Channel1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.


Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human kidney.

CLIC1



See Immunoblot Data and Customer Feedback for more information

Antibody used in WB on Mouse renal epith...



Antibody used in WB on Human Jurkat 1.25...