

Product Datasheet

SCNN1B Antibody (orb1243218)

Description

SCNN1B Antibody

Species/Host

Rabbit

Reactivity

Canine, Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human SCNN1B.

Antibody used in WB on Human Jurkat 0.12...

Target

SCNN1B

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SCNN1B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in WB on Human A549, H441.

Note

For research use only

Application notes

SCNN1B antibody can be used for detection of SCNN1B by ELISA at 1:312500. SCNN1B antibody can be used for detection of SCNN1B by western blot at 0.125 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

70 kDa

Uniprot ID

[P51168](#)

NCBI

[NP_000327](#)

Dilution Range

SCNN1B antibody can be used for detection of SCNN1B by ELISA at 1:312500. SCNN1B antibody can be used for detection of SCNN1B by western blot at 0.125 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90kDa

80kDa

70kDa

60kDa

50kDa

40kDa

30kDa

20kDa

10kDa

5kDa

2kDa

1kDa

500Da

250Da

125Da

62Da

31Da

15Da

7Da

3Da

1Da

500Da

250Da

125Da

62Da

31Da

15Da

7Da

3Da

1Da

500Da

250Da

125Da

62Da

31Da

15Da

7Da

3Da

1Da

500Da

250Da

125Da

62Da

31Da

15Da

7Da

3Da

1Da

500Da

250Da

125Da

62Da

31Da

15Da

7Da

3Da

1Da

500Da

250Da

125Da

62Da

31Da

15Da

7Da