

Product Datasheet

CLCNKB Antibody (orb1243219)

Description

CLCNKB Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
kidney at 0...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CLCNKB.

Target

CLCNKB

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CLCNKB antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CLCNKB antibody can be used for detection of CLCNKB by ELISA at 1:312500. CLCNKB antibody can be used for detection of CLCNKB by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

76 kDa

Uniprot ID

[P51801](#)

NCBI

[NP_000076](#)

Dilution Range

CLCNKB antibody can be used for detection of CLCNKB by ELISA at 1:312500. CLCNKB antibody can be used for detection of CLCNKB by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.