

Product Datasheet

CLCNKB Antibody (orb1243220)

Description

CLCNKB Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CLCNKB.
Target	CLCNKB
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store CLCNKB antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	CLCNKB antibody can be used for detection of CLCNKB by ELISA at 1:312500. CLCNKB antibody can be used for detection of CLCNKB by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	50 kDa, 75 kDa
Uniprot ID	Q8WW53
NCBI	AAH20873
Dilution Range	CLCNKB antibody can be used for detection of CLCNKB by ELISA at 1:312500. CLCNKB antibody can be used for detection of CLCNKB by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Jurkat 0.2-...