

Product Datasheet

ZMYM6 Antibody (orb1243234)

Description

ZMYM6 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZNF258.

Target

ZMYM6

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZNF258 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ZNF258 antibody can be used for detection of ZNF258 by ELISA at 1:62500. ZNF258 antibody can be used for detection of ZNF258 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

79 kDa

Uniprot ID

[O95789](#)

NCBI


[NP_009098](#)

Dilution Range


ZNF258 antibody can be used for detection of ZNF258 by ELISA at 1:62500. ZNF258 antibody can be used for detection of ZNF258 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in WB on Human Jurkat 0.2-...