

Product Datasheet

CREM Antibody (orb1243247)

Description

CREM Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CREM.

Target

CREM

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CREM antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CREM antibody can be used for detection of CREM by ELISA at 1:1562500. CREM antibody can be used for detection of CREM by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

13 kDa

Uniprot ID

[O75519](#)

NCBI


[NP_874386](#)

Dilution Range

CREM antibody can be used for detection of CREM by ELISA at 1:1562500. CREM antibody can be used for detection of CREM by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human Jurkat 2.5 ...