

Product Datasheet

KLHL13 Antibody (orb1243280)

Description

KLHL13 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KLHL13.

Target

KLHL13

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store KLHL13 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

KLHL13 antibody can be used for detection of KLHL13 by ELISA at 1:62500. KLHL13 antibody can be used for detection of KLHL13 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

74 kDa

Uniprot ID

[Q9P2N7](#)

NCBI


[NP_277030](#)

Dilution Range


KLHL13 antibody can be used for detection of KLHL13 by ELISA at 1:62500. KLHL13 antibody can be used for detection of KLHL13 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Intestine.



Antibody used in WB on Human Jurkat 0.2-...