


Product Datasheet

DMTF1 Antibody (orb1243296)


Description

DMTF1 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human DMTF1.
Target	DMTF1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store DMTF1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	DMTF1 antibody can be used for detection of DMTF1 by ELISA at 1:62500. DMTF1 antibody can be used for detection of DMTF1 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	54 kDa, 84 kDa
Uniprot ID	Q9Y222
NCBI	EAW76962
Dilution Range	DMTF1 antibody can be used for detection of DMTF1 by ELISA at 1:62500. DMTF1 antibody can be used for detection of DMTF1 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Intestine.



Antibody used in WB on Human HepG2 at 0....