

Product Datasheet

ZNHIT6 Antibody (orb1243303)

Description

ZNHIT6 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human FLJ20729.

Target

ZNHIT6

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store FLJ20729 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

FLJ20729 antibody can be used for detection of FLJ20729 by ELISA at 1:312500. FLJ20729 antibody can be used for detection of FLJ20729 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

54 kDa

Uniprot ID

[Q9NWN0](#)

NCBI


[NP_060423](#)

Dilution Range


FLJ20729 antibody can be used for detection of FLJ20729 by ELISA at 1:312500. FLJ20729 antibody can be used for detection of FLJ20729 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Intestine.



Antibody used in WB on Human Jurkat 1.25...