

Product Datasheet

TBC1D10B Antibody (orb1243320)

Description		
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TBC1D10B.	Antibody used in WB on Human HepG2 at 1:....
Target	TBC1D10B	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store TBC1D10B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	TBC1D10B antibody can be used for detection of TBC1D10B by ELISA at 1:12500. TBC1D10B antibody can be used for detection of TBC1D10B by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	61 kDa	
Uniprot ID	Q9H8Z2	
NCBI	NP_056342	
Dilution Range	TBC1D10B antibody can be used for detection of TBC1D10B by ELISA at 1:12500. TBC1D10B antibody can be used for detection of TBC1D10B by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	