

Product Datasheet

COLQ Antibody (orb1243370)

Description	COLQ Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human COLQ.	
Target	COLQ	Antibody used in WB on Human HepG2 at 0.... 
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	Antibody used in WB on Hum. Adult Placen... 
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store COLQ antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	COLQ antibody can be used for detection of COLQ by ELISA at 1:62500. COLQ antibody can be used for detection of COLQ by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	34 kDa	
Uniprot ID	Q9Y215	Antibody used in WB on Hum. Adult Placen... 
NCBI	NP_536805	
Dilution Range	COLQ antibody can be used for detection of COLQ by ELISA at 1:62500. COLQ antibody can be used for detection of COLQ by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	