

Product Datasheet

TMEM260 Antibody (orb1243375)

Description

TMEM260 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human C14orf101.

Target

TMEM260

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store C14orf101 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

C14orf101 antibody can be used for detection of C14orf101 by ELISA at 1:12500. C14orf101 antibody can be used for detection of C14orf101 by western blot at 0.125 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

79 kDa

Uniprot ID

[Q9NX78](#)

NCBI

[NP_060269](#)

Dilution Range

C14orf101 antibody can be used for detection of C14orf101 by ELISA at 1:12500. C14orf101 antibody can be used for detection of C14orf101 by western blot at 0.125 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human HepG2 at 0....