

Product Datasheet


CCT4 Antibody (orb1243395)

Description

CCT4 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CCT4.
Target	CCT4
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store CCT4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	CCT4 antibody can be used for detection of CCT4 by ELISA at 1:62500. CCT4 antibody can be used for detection of CCT4 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	58 kDa
Uniprot ID	Q53QP9
NCBI	NP_006421
Dilution Range	CCT4 antibody can be used for detection of CCT4 by ELISA at 1:62500. CCT4 antibody can be used for detection of CCT4 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

CCT4




White: CCT4


See H2C.1 Data and Customer Feedback for more information

Antibody used in IHC on Rat hippocampal

...



Antibody used in WB on Human HepG2 at 0....



Antibody used in WB on Human HEK 293T at...