


Product Datasheet

PREB Antibody (orb1243425)

Description	PREB Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PREB.	Antibody used in IHC on Human Heart.
Target	PREB	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PREB antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in WB on Human Raji at 0.2...
Note	For research use only	
Application notes	PREB antibody can be used for detection of PREB by ELISA at 1:312500. PREB antibody can be used for detection of PREB by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	45 kDa	
Uniprot ID	Q9HCU5	
NCBI	NP_037520	
Dilution Range	PREB antibody can be used for detection of PREB by ELISA at 1:312500. PREB antibody can be used for detection of PREB by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	