

Product Datasheet

UBE2N Antibody (orb1243436)

Description

UBE2N Antibody

UBE2N

Species/Host

Rabbit

Reactivity

Canine, Drosophila, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human UBE2N.

Target

UBE2N

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store UBE2N antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

UBE2N antibody can be used for detection of UBE2N by ELISA at 1:312500. UBE2N antibody can be used for detection of UBE2N by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

17 kDa

Uniprot ID

P61088

NCBI


NP_003339

Dilution Range

UBE2N antibody can be used for detection of UBE2N by ELISA at 1:312500. UBE2N antibody can be used for detection of UBE2N by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.


Expiration Date

12 months from date of receipt.



See Immunoblot 3 Data and customer Feedback for more information

Antibody used in WB on r-UBE2N at 1:500 ...



See IHC 1 Data and Customer Feedback for more information

Antibody used in IHC on Human brainstem ...



Antibody used in WB on Human HepG2 cells...