


Product Datasheet

PARP6 Antibody (orb1243437)

Description

PARP6 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PARP6.
Target	PARP6
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PARP6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	PARP6 antibody can be used for detection of PARP6 by ELISA at 1:1562500. PARP6 antibody can be used for detection of PARP6 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	59 kDa
Uniprot ID	Q2NL67
NCBI	NP_064599
Dilution Range	PARP6 antibody can be used for detection of PARP6 by ELISA at 1:1562500. PARP6 antibody can be used for detection of PARP6 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Intestine.



Antibody used in WB on Human HepG2 at 0....