

Product Datasheet

BHLHE40 Antibody (orb1243447)

Description

BHLHE40 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human BHLHB2.

Target

BHLHE40

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store BHLHB2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

BHLHB2 antibody can be used for detection of BHLHB2 by ELISA at 1:62500. BHLHB2 antibody can be used for detection of BHLHB2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

45 kDa

Uniprot ID

[O14503](#)

NCBI

[NP_003661](#)

Dilution Range


BHLHB2 antibody can be used for detection of BHLHB2 by ELISA at 1:62500. BHLHB2 antibody can be used for detection of BHLHB2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in WB on Transfected 293T ...



Antibody used in WB on Hum. Fetal Lung a...



Antibody used in WB on Hum. Adult Placen...