



Product Datasheet

CNP Antibody (orb1243462)

Description	CNP Antibody	CNP
Species/Host	Rabbit	
Reactivity	Human	Antibody used in IP on Mouse Brain.
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CNP.	
Target	CNP	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	Antibody used in WB on Human NT-2, Rat b...
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store CNP antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	CNP antibody can be used for detection of CNP by ELISA at 1:62500. CNP antibody can be used for detection of CNP by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	45 kDa	Antibody used in WB on Human Jurkat 1.25...
Uniprot ID	P09543	
NCBI	NP_149124	
Dilution Range	CNP antibody can be used for detection of CNP by ELISA at 1:62500. CNP antibody can be used for detection of CNP by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	