

Product Datasheet

SMPDL3B Antibody (orb1243465)

Description

SMPDL3B Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SMPDL3B.

Target

SMPDL3B

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SMPDL3B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SMPDL3B antibody can be used for detection of SMPDL3B by ELISA at 1:312500. SMPDL3B antibody can be used for detection of SMPDL3B by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

52 kDa, 42 kDa

Uniprot ID

[Q92485](#)

NCBI


[NP_055289](#)

Dilution Range

SMPDL3B antibody can be used for detection of SMPDL3B by ELISA at 1:312500. SMPDL3B antibody can be used for detection of SMPDL3B by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Intestine.



Antibody used in WB on Human Raji at 0.2...